




NESSY

NMR Relaxation Dispersion Spectroscopy Analysis Software for Windows, Mac and Linux



Manual

Functions are grouped according to the menu entries:

Menu	Description of function
------	-------------------------

File:

New	Empties all settings and variables to start a new NESSY project.
Open	Opens a saved NESSY project.
Save	Saves a NESSY project.
Save as	Saves a NESSY project under new name.
Opened Projects	List of NESSY projects that were loaded previously.
Quit	Closes NESSY. Note that running calculations will be aborted.

Import:


Load Structure File

Reads the sequence of a PDB file. Numbering will be identical to numbering of PDB structure.

It is recommended to read in the sequence from a PDB file, as it is the most convenient way to have it synchronized with your data and results. Note that only colour and width coded structures can be viewed if a PDB file was selected/imported.

Load FASTA sequence

Opens up a dialog to add protein sequence in FASTA format.



The sequence can be retrieved from the internet by indicating UniProtKB/Swiss-Prot code (www.expasy.ch) or by writing/pasting it in the panel. If the sequence is not starting with residue no. 1 of the current project, this can be set.

By pushing the clear button, the loaded sequence will be deleted from the Data tab.

Import Data from Single File

Use this function to import a peak list to your NESSY session. First, you have to load the peak list (Datafile). Then you can specify in which data set row the data should be imported. NESSY is able to import (almost) all peak lists. Therefore, the column number of the assignment (column where residue number is) and the intensity has to be specified. Note that numbering starts with the first real entry. Finally, you can choose between different peak file formats. Sparky, Xeasy and NMR View peak file formats can be selected. In addition, you have the opportunity to indicate whether your peak list is a space/tabulator or semicolon separated text file.

Example (Sparky peak list):

Assignment	w1	w2	Data		Sequence	Data 1
D3N-H	122.284	8.519	6268538	3	Asp	6268538
C4N-H	116.715	7.968	98160120	4	Cys	98160120
I5N-H	124.403	8.893	39739664	5	Ile	39739664
K7N-H	116.092	7.58	39971648	7	Lys	39971648
W8N-H	115.785	9.255	7588816	8	Trp	7588816
K9N-H	119.422	7.735	20583350	9	Lys	20583350
1	2	3	4			

Column Numbering

Import Data from Multiple Files


Use this function to import multiple peak lists at once with the same format. Set up is similar than Import Data from Single File, except you can select multiple files and you specify the dataset from where the data will be appended.

Import VD List

This function allows you to import the vd list of your (Bruker) CPMG experiment. CPMG frequencies are listed in the vd list and will be imported directly in your NESSY session.

Import BRUKER Protein Dynamic Center Project

This function is used to import projects started with BRUKER's Protein Dynamics Center (PDC).



Use PDC to export the project as text file (see PDC manual) and select this file in the dialog shown above. If you are using multiple experiments (for global fit), you can specify in which project the data has to be loaded.


Data:

Correction of Residue Numbering	Shifts your sequence up or downwards.
Correction of Data Assignment	Shifts you data up or downwards.
Erase Data Set	Function to delete single data sets.

Settings:

Select residues to include in calculation


Opens a dialog where users can select/deselect residues, which are included in further analysis.



Residues previously loaded into the project (by FASTA or pdb) can be selected. Entire protein, only few or single residues can be selected.

Select Models

opens a dialog to enable/disable individual models.




By default, models 1 - 3 are active. Each model, except of model 1, can be switched off or on. If multiple experiments are set up, data can be globally fit. To do so, Global Fit model has to be activated. Now, data is only globally fitted to selected models. If multiple experiments are set up, but Global Fit is not activated, experiments are individually fitted.

Note: Selection of models is saved in NESSY projects, but Global Fit has to be activated after loading the experiments!

NESSY Settings

You can set properties for your data analysis:



- Model Selection (Default is AICc):
 - Akaike Information criteria with a second order for small sample size (AICc)
 - Akaike Information criteria (AIC)
 - F-Test
 - Chi2
- Number of Monte Carlo Simulations (default is 500)
- Maximum number of data sets
- Minimal difference between R2eff between lowest or highest CPMG frequency in Hz.

This setting can be used to predict chemical exchange and to speed up calculation. Default difference is 2 Hz.


- Calculate models 2-5 lets users choose if selected models are calculated for each residue ('for entire dataset') or only if chemical exchange is expected ('only if exchange is expected').
- Range of Y axis for R2eff plots (including fits)
- Width of frame for generated plots
- Type of vector graphics (plots), which are generated during calculation. Supported formats are: SVG, PS, EPS, PDF (default is SVG)

NESSY uses AICc model selection as default, as normal distribution of data can not be assumed (as necessary for F-Tests). Nevertheless, settings can be changed and defaults can be resored at any time of and during your calculation. Note that if you change settings during the calculation, your data might not be uniformly analyzed.

n-States:

Fast Exchange


Opens up dialog to perform fit to 2-9 states models in fast exchange.



Residues to analyze and limit for number of states can be specified.

Slow Exchange

Opens up dialog to perform fit to 2-9 states models in slow exchange.



Residues to analyze and limit for number of states can be specified.

Extras:

PEAKY

Automatic peak picking and assignment tool (*not implemented yet !*).

Back calculate intensities

Calculates and normalizes dispersion data of independent experiment to your session to implement additional data sets.

This function can be used for:


- Adding data points of independent performed experiments
- Combine data recorded at different spectrometer frequencies

The main function is to make a baseline correction of independent experiment so data can be fitted together to dispersion models. This functions generates artificial Sparky-formatted peak lists that can be imported into your current session.

NESSY calculates R2eff of your additional datasets and back calculates intensities (normalized according to reference spectra of your NESSY session). This allows to import and combine data of different experiments.

Create Synthetic data

This function can be used to create synthetic data for models 2-5.



CPMG frequencies can be specified (separated by comma) as well as the constant time for the CPMG pulse train. Synthetic peak lists (in Sparky format) will be generated and saved in

selected folder for specified model. The error, which will be introduced, can be set using the slider.

Calculate free Energy dG

Free energy can be calculated out of relaxation dispersion experiments.

By inserting dynamics parameters calculated by NESSY, free energy (dG) can be calculated. If model 2 (fast exchange) was selected, populations can be calculated if dw (shift differences) are known (e.g. for binding experiments).

Calculated dG values can be directly passed to van't Hoff analysis (see below).

Formula:

$$\Delta G = RT \ln \left(\frac{p_a}{p_b} \right)$$

van't Hoff Analysis

If relaxation dispersion experiments are performed at different temperatures, van't Hoff analysis can be used to calculate Enthalpy and Entropy.

The linear and non-linear van't Hoff models are included for calculation. Either dG (free energy) can be added to data grid (manually or directly after calculation; see above) or K (= $k(b-a) / k(a-b)$) together with the temperature in Kelvin.

Data does not have to be obtained by NESSY calculation. This function can be used to perform van't Hoff analysis for different experiments.


Export/Save log

Execution protocol can be saved in a text file to track calculation.

Plotting:

Create 2D plot

This function can be used to create custom made 2D plots.



Either a NESSY csv file or any other text file (comma or semicolon separated columns) can be loaded into dialog. For each data set, a figure legend can be indicated. In addition, users have the choice between 5 different plot styles and can select the color for the data set. The column number of x-, y- and error-values are specified per data set (column numbering starts with 1 for the first column; white space does not count).


In addition, labels of x and y axis as well as the range can be set.

A file root for the plots has to be given. Two files are created: a png file, which will be added to the results tab in the NESSY main windows, and a vector graphics plot (depending on the selection in NESSY settings).

Note that NESSY plot generator can be used to plot data not obtained by NESSY. Any data in text files can be inserted (comma or semicolon separated data). Location of important data (x, y or error) can be selected.

Create 3D plot


This function can be used to create custom made 3D plots.



3D plot generator works identical as the 2D plot generator. NESSY or any other text file can be used as source. The only difference is that the 3rd dimension can be specified using the 'Plane no.' text field. Multiple data sets can have the same plane number. This is useful to plot a scatter plot with data points and a linear plot, which is the regression, in the same plane.

Create Multiple Alignments of Plots

This function can be used to create a SVG graphic of multiple aligned plots.




The number of columns can be selected. NESSY calculates the corresponding numbers of rows.

If a NESSY project is loaded or data was analyzed, plots created during analysis will be added to the list in the dialog shown above. Alternatively, SVG files can be drag&dropped into the plots field.

A SVG graphics of the plots (or graphics selected) will be created.

Note that plots are sorted by increasing residue number!



This function is very useful to generate an ensemble of individual fits for each residue to present in supplementary materials. Time to generate the plot above was less than 1 minute!

Help:

Manual

This manual will appear.

Tutorial

Follow the tutorial to learn how to use NESSY's basic functions.

GNU GPL License

Read the GPL license of NESSY.

NESSY mailing list

Join the NESSY mailing list on <http://gna.org/mail/?group=nessy>

Contact NESSY

Write an email to the NESSY mailing list (using default mail program).

Update

If you installed the SVN version of NESSY, you can synchronize NESSY with further developments.

For further questions or support please subscribe at our [mailing list](#).

(C) Michael Bieri 2010